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(1) Two-step glutaraldehyde method (A 產品)

1. 10mg horseradish peroxidase (RZ 3.0) is dissolved in 0.2ml 0.1M PBS pH 6.8 containing 1.25% glutaraldehyde. Mixture is left overnight at room temperature.
2. The mixture is dialyzed against normal saline or passed down a Sephadex G25 column to remove free glutaraldehyde. Solution is made up to 1.0ml in normal saline.
3. The globulin solution of the anti-serum is dialyzed against PBS at +4°C and is adjusted to 5mg/ml in normal saline.
4. 1.0ml of the globulin solution is mixed with 1.0ml of the activated HRP solution and 0.1ml 1M carbonate bicarbonate buffer pH 9.5 is added and the mixture is allowed to stand at +4°C for 24 hours.
5. 0.1ml 0.2M lysine solution is added and mixture is kept at room temperature for 2 hours.
6. The mixture is dialyzed against PBS pH 7.2 overnight.
7. The enzyme-labelled globulin is precipitated by addition of an

equal volume of saturated ammonium sulphate solution. Precipitate is washed twice in half saturated ammonium sulphate, and is suspended in 1.0ml PBS.

8. The conjugate is dialyzed against PBS extensively (or passed down a Sephadex G25 column).
9. Conjugate is centrifuged at 10,000G for 30 minutes and sediment is discarded. Bovine serum albumin is added to 1%.
10. Conjugate is filtrated through a millipore filter (0.22µm). It can be stored at -20°C or at +4°C if an equal amount of glycerol is added.

(2) Periodate method (Wilson & Nakane) (B 產品)

Suitable for peroxidase.

1. Dissolve 4mg horse-radish peroxidase (RZ 3.0) in 1.0ml H₂O. Add 0.2ml freshly made 0.1M NaIO₄. Stir for 20 minutes at room temperature.
2. Dialyze against 1mM sodium acetate buffer pH4.4 overnight at +4°C.
3. Add 20ul 0.2M carbonate buffer pH 9.5 and at once add 8mg immunoglobulin dissolved in 1.0ml 0.01M carbonate buffer pH 9.5. Stir 2 hours at room temperature.
4. Add 0.1ml fresh sodium borohydride solution (4mg/ml in H₂O). Allow to stand 2 hours at +4°C.
5. Dialyzed against PBS, then repeated steps 9 and 10 in method (1).

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3. Avrameas, S., and B. Guibert 1971 Biologically active waterinsoluble protein polymers. Their use for isolation of specific interaction proteins. *Biochemistry* 53:603-629.
- 4 Lanner, M., R. Berguist, J. Carlsson, G. Huldt 1978. Purification of Enzyme-labelled conjugation by affinity chromatography. p.237-241, in "Affinity Chromatography" ed. Hoffmann-Ostenhof et al. Pergamon Press. Oxford.
5. Voller, A., D.E. Bidwell, and A. Bartlett 1979. The enzyme linked immunosorbent assay (ELISA). p.39-40. Nuffield laboratories of Comparative Medicine. The Zoological Society of London.

THE DEVELOPMENT OF ENZYME CONJUGATION AND ITS TECHNIC IMPROVEMENT

Fei, C.Y., T.S. Huang, P.H. Mar*

Taiwan Provincial Research Institute for Animal Health.

Rabbit anti-swine IgG antibody was conjugated with horse raddish peroxidase by the glutaldehyde and Nakane methods to evaluate the effect of the two conjugation methods.

Results showed the product of the glutaldehyde method could be further purified by ConA-sepharose to separate the free antibody and the conjugated antibody. The product of the Nakane method could not be further purified by the ConA-sepharose because the sugar of the enzyme was oxygenated in the conjugation process. The molarity ratio of enzyme/antibody of the glutaldehyde and Nakane methods were 1.5 and 0.9 respectively.

* National Taipei College of Nursing