STELLA KO
MTTW, 2022

USP643 and TOC Theory

Total Organic Carbon

METTLER TOLEDO

■ USP

- Relevant Sections of the USP
- 〈643〉 TOC - Calibration & Performance
- 〈643〉 TOC – SST
- TOC Maintenance

■ Theory

- Why Measure TOC?
- TOC Basic Theory
- Principals of Operation



■ Application

- UPW Make-up & Distribution System
- Pharma Application

■ Instrument


- Introduction
- Workshop

- The USP is the U.S. Pharmacopeia, the latest version is USP 32
- It is a private, not-for-profit, *non-governmental* organization.
- There are no political or governmental links.
- **It is dedicated to the development of quality standards for the benefit of the public health.**
- To improve the health of people around the world through public standards and related programs that help ensure the **Quality, Safety,** and **Benefit** of medicines and foods.

The USP establishes standards for the benefit of public health... they are legally enforced by the FDA.

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Bulk Water for Injection(s) - Today

Attribute ¹	USP	EP	JP
Production Method	Distillation or suitable process	Distillation	Distillation, RO with UF, from Purified Water
Source Water	US, EU, Japan, WHO drinking water	Human consumption	JP water specification
Total Aerobic (cfu/100 mL) ²	-	10	-
Conductivity (μ S/cm at 25° C) ³	1.3 (3 stage)	1.3 (3 stage)	2.1 offline
TOC (mg/L)	0.5	0.5	0.5 (0.3 for control)
Bacterial Endotoxins (EU/mL)	0.25	0.25	0.25
Nitrates (ppm)		0.2	Effective JP 16
Heavy Metals (ppm)			
Acidity/Alkalinity			
Chloride			
Sulfate			
Nitrite			
Ammonium (mg/mL)			
Oxidizable Substances (/100 mL)			
Residue on Evaporation (mg)			

Note 1: All tests are maximum, unless otherwise stated.

Note 2: Microbiological testing is considered to be harmonized, with the exception noted that the EP test is written into the Production section, and the USP test is contained in a non-compendial general information chapter

Note 3: Limits are temperature dependent

Note 4: Not required effective Jan 1, 2009.

Attribute ¹	USP	EP	JP
Production Method	Suitable process	Suitable process	Distillation, ion-exchange, UF, or combination
Source Water	US, EU, Japan, WHO drinking water	Human consumption	JP water specification
Total Aerobic (cfu/mL) ²	100	100	100
Conductivity (μ S/cm at 25° C) ³	1.3 (3 stage)	5.1 (1 stage)	2.1 offline
TOC (mg/L)	0.5	0.5 (optional)	0.5
Bacterial Endotoxins (EU/mL)			
Nitrates (ppm)		0.2	Effective JP 16
Heavy Metals (ppm)			
Acidity/Alkalinity			
Chloride			
Sulfate			
Nitrite			
Ammonium (mg/mL)			
Oxidizable Substances (/100 mL)		<0.1 mL ⁴ 0.02 KMnO ₄	
Residue on Evaporation (mg)			

Note 1: All tests are maximum, unless otherwise stated.

Note 2: Microbiological testing is considered to be harmonized, with the exception noted that the EP test is written into the Production section, and the USP test is contained in a non-compendial general information chapter

Note 3: Limits are temperature dependent



Note 4: Alternative to TOC

Note 5: Not required effective Jan 1, 2009 if WFI conductivity requirements are met


- Limit of Detection of 0.050 mg C/L (50 ppb TOC)
- Must distinguish inorganic carbon, i.e., CO_2 , HCO_3^-
- Must meet System Suitability testing periodically
- Calibrate according to the manufacturer's recommendations

- **Question:** What is System Suitability Testing?
- **Answer:** The process of challenging a TOC analyzer to convert sucrose and p-benzoquinone to another detectable species (CO_2) equally.

It is required by USP <643> and EP 2.2.44

■ 1,4-benzoquinone chemical structure

■ Sucrose molecule chemical structure

- Measure TOC of 0.50 mg C/L (500 ppb C as sucrose), R_s .
- Measure TOC of 0.50 mg C/L (500 ppb C as *p*-benzoquinone), R_{ss} .
- Measure TOC of water used to prepare these solutions, R_w . Not to exceed 100 ppb.
- Response shall be between 85 and 115%!

$$\text{Response} = 100 \times \left(\frac{R_{ss} - R_w}{R_s - R_w} \right)$$

- R_s , Sucrose ~516 ppb
- R_{ss} , p-benzoquinone ~499
- R_w , Reagent Water ~15-20 ppb

$$\text{Response} = 100 \times \left(\frac{499 - 15}{516 - 15} \right) = 96.7\%$$

Standard Service Products -

After installation, service required every 6 months

*For Pharmaceutical end-users only

Basic Preventative Maintenance	UV Lamp and Filter Replacement TOC Calibration
Full Preventative Maintenance	UV Lamp and Filter Replacement Flow, Temperature, Conductivity and TOC Calibration

- **Pharma – 500 ppb** Residual compounds remain in the manufacturing process, the potential for side effects from toxic levels of contaminants increases
- **MicroE – 2-3 ppb** Oxidation film has poor insulation withstand voltage and reduces Wafer's life
- **Power – 20ppb** The high temperature and high pressure environment of the generator oxidizes organic carbon molecules into carbonic acid, which in turn erodes the turbine

Water	Typical TOC (Range)	Typical TOC present as Particulates
Bog	33 ppm (10 to 60)	4 ppm
Marsh	17 ppm (10 to 60)	3 ppm
Eutropic lake	12 ppm	3 ppm
Oligotrophic lake	2.2 ppm	0.2 ppm
River	7.0 ppm (1 to 10)	2.5 ppm
Precipitation	1.1 ppm	0.1 ppm
Ground Water	700 ppb	
Sea Water	500 ppb	50 ppb
Waste Water	up to 1000 ppm	
Process Waters	very wide range	
Drinking Water	100 ppb to 10 ppm	
Purified Water	1 ppb to 500 ppb	
Ultrapure Water	0.1 ppb to 10 ppb	

Table Showing Total Organic Carbon levels in various types of water

- Leachates and humic acids and other organic sources
 - Large molecules
 - Complex mixtures
 - Some man-made waste... pesticides?

leaves

most
probable


leachate (*n.*) A product or solution formed by leaching, especially a solution containing contaminants picked up through the leaching of soil. Pesticides or fertilizers can leach through soil.

Humic acid is one of the major components of humic substances which are dark brown and major constituents of soil organic matter humus that contributes to soil chemical and physical quality and are also precursors of some fossil fuels. They can also be breakdown products from **leaves**.

- Degradation of the water system
 - RO membranes
 - Filters
 - Ion-exchange resins
 - High Purity Water Pumps
- Formation within water system
 - Biofilm?


particulates

biofilm


A **biofilm** is a complex aggregation of microorganisms marked by the excretion of a protective and adhesive matrix.

- TOC is about the ability to break bonds and to convert "undetectable organic" molecules into a chemical form that is detectable
 - Convert to methane in a GC/FID (Flame Ionization Detector)
 - Convert to CO₂ and detect by infrared
 - Convert to CO₂ and detect by conductivity
- This is done by BREAKING BONDS and OXIDIZING
- UV oxidation breaks the carbon-carbon bonds through the generation of hydroxyl radicals to produce carbon dioxide and water. It is the dissolved carbon dioxide that we then measure.


■ **Sucrose molecule
chemical structure**


■ **1,4-benzoquinone
chemical structure**


- Conductivity measurement in UPW is the measurement of total ionic (conductive) impurity.
 - Dozens of ions in water
 - No specificity
- TOC measurement in UPW is the measurement of total organic (usually non-conductive) impurity.
 - Hundreds/thousands of organics in water
 - No specificity
- TOC is a conductivity based measurement. Thornton is a world leader in this field

On-Line 6000TOC Measurements by UV Oxidation/Conductivity


water)

$$\text{Outlet (TOC+TIC)} - \text{Inlet (TIC)} = \text{TOC}$$


Ultra violet light at 185 nm breakdowns non-ionic organic species which are then measured by conductivity sensors

The Same continuous flow and rapid response.


Control valve automatically maintains flow rate!

- 1964 Dr. Richard Thornton, an MIT professor, founded Thornton Inc. in USA
2001 Acquired by Mettler-Toledo and integrated into Process Analytics Division

Widely considered as the leader in the production of high quality resistivity / conductivity sensors and the Thornton/Light Curve
Long history of community marketing successful at establishing brand recognition at key international organizations such as USP (A. Bevilacqua), EP, JP, FDA, SEMI & PV (J. Cannon) and ASTM (D. Gray)
Strong working relationships at system fabricators

The audit-ready solution for all TOC applications

Water Purification

- Monitor TOC levels of systems producing:
 - Water for Injection (WFI)
 - Purified Water (PW)

Final Rinse

- Ensure **final rinse** completion with TOC after Clean in Place (CIP) applications

Point-of-use

- TOC sensor located after RO purification, storage and distribution at water point-of-use

What can happen in 6 minutes?



- A lot of good (and bad) measurement data can be missed in 6 minutes - if your data is limited to snapshots of information
- Unlike at-line batch analyzers, the 6000TOCi provide **continuous, real-time data** to show exactly when an excursion starts and ends.


Watch the video*
"What can happen in 6 minutes?"

6000TOCi Sensor - Continuous, real-time measurement

At-line Batch Analyzer - Batch measurement, 6+ minutes

VS.

*Internet connection required to watch video

Verifiable sensor performance with Predictive Diagnostics

- The 6000TOCi with the M800 transmitter offers a unique set of diagnostics to verify during an audit that your TOC analyzer was operating appropriately during its use.
- Ensures your TOC system is always measuring effectively and efficiently.

Lamp Dynamic Lifetime Indicator

Predictive diagnostic tool that ensures the UV oxidation lamp is providing oxidation according to specification.

TTM – Time to Maintenance

Indicates when services are required, including Calibration, SST and ballast.

Sensor Alarms & Excursion Counters


Clearly alerts you if there is an identifiable problem with your TOC measurements.

Advanced diagnostic tool for simple maintenance

- The lamp DLI predictive diagnostic tool ensures accurate status of your lamp life hours remaining to minimize unplanned maintenance
- Based on known lamp life characteristics the lamp DLI algorithm projects the lamp lifetime based on power cycling of the TOC lamp
- All UV lamps degrade over time. Much like how cycling a light switch on and off degrades the life of your household lights, so does power cycling of your TOC lamp
 - Multiple lamp power cycles per hour rapidly degrades a lamp
 - Few lamp power cycles per day has very little impact
 - Lamp DLI can distinguish this difference

Solid design based on proven technology is more robust with improved oxidation performance characteristics

Stainless steel fittings

Advanced optical modeling enhancements improve repeatability, reproducibility and stability

Fast (< 1 minute response time)

New, more robust coil with oxidation chamber view port

Durable SS tube fittings

Viewing port allows quick inspection for contamination or rouge

A robust product meeting process application requirements

Rugged Enclosure

- Powder coated metal back enclosure
- Molded plastic front enclosure
- Easy to clean
- Consistent PRO Branding

USB

- Direct printer capabilities
- Data logging

Service Access

- Full access through front cover
- Side access for UV lamp replacement

RFID Technology

- Semi-automated calibration and system suitability test
- Radio frequency identification (RFID) enabled for error proof solution data transfer

LED Indicators

- At-a-glance status

Isolated Connections

- Left side electrical connections
- Right side water connections

A functional and effective design for dependable measurement

Robust Enclosure

- Powder coated metal
- Gasket lined

Separated Electrical & Wet Components

- Elevated PCB & electrical connectors

Removable Service Panels


- PCB board
- Back panel

Easy Access Side panel

- UV lamp access
- USB printer capable
- USB data logger

Universal Ballast

- 100-240 VAC

Centralized Signal & Control Connectors

Efficient Sample Flow Path


- Compression fittings
- 316SS tubing
- Shortened flow path

Conductivity Sensor Enhancement

- Optimized flow velocity for more efficient measurement

Oxidation Chamber Enhancement

- Quartz coil optimization
- Quartz chamber view port

Proven technology - NEW features and value for improved TOC performance!

Consistent Measurements with Stability

Audit-ready continuous, Real-time TOC Analysis

Enhanced Data Assurance for Global Compliance

Automatic solution recognition and data transfer minimizes errors and saves time

RFID identification requires MT solutions

RFID tag ensures solutions are installed in the correct order

Smaller 250 ml bottles are more convenient and easy to handle

RFID automatically transmits unique solutions data directly to M800

450TOC

- Portable analyzer
- USB outputs
- Full regulatory compliance

4000TOCe with M300 G1

- Basic TOC measurement
- Full regulatory compliance
- Transmitter Configuration Tool (TCT)

6000TOCi with M800

- Premium sensor and powerful touch-screen M800 transmitter and TCT
- TOC performance enhancements and stability
- USB data logging and printer capabilities
- New ISM diagnostics

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The screenshot shows the Mettler-Toledo website. At the top, there is a navigation bar with links for Home, Products & Solutions, Services & Support, and About Us. Below the navigation bar, the main content area features a portrait of Jim Cannon, a Pharmaceutical Waters Industry Expert. A quote asks, "How often must we calibrate conductivity sensors?" and explains that in the pharmaceutical industry, one year is the most accepted time frame. Below the quote, there is a section titled "Product Solutions" with links to UVCond Conductivity/Resistivity, UVCond Calibrators, Total Organic Carbon, Dissolved Ozone and Oxygen, Seawater Analysis, Silica Analysis, and Transmitters. There are also three columns of questions and answers related to regulations, applications, and calibration.

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© 03/14 Mettler-Toledo Thornton, Inc.
Mettler-Toledo Pac Rim AG, Taiwan Branch